

# Cleaning and Disinfecting Guide

## Disinfectant List for DR Systems

- The following disinfectants can be used to disinfect appearance & cover of digital X-ray systems  
: GC85A, GM85, GC80, GC70

Name	Manufacture	Active Ingredient	CBC*	EPA*
Accel PREvention Wipes	Diversey Canada, Inc.	Hydrogen peroxide		
ALI-FLEX RTU	Lalema	Sodium hypochlorite		
ANIOSURF Premium NPC	ANIOS	Quaternary Ammonium(N-Alkyl)		
ANIOS社 surfa'safe	ANIOS	Quaternary Ammonium(N-Alkyl)		
Cleanisept Wipes MAXI	Dr.Schumacher GmbH	Quaternary Ammonium(N-Alkyl)		
CLOROX Bleach Germicidal	The Chlorox Company	Sodium hypochlorite	○	○
Dettol	Reckitt Benckiser	Quaternary Ammonium(N-Alkyl)		
Ecolab Actichlor plus	Ecolab	Troclosene sodium		
ED Wipes	MH healthcare	Quaternary Ammonium(N-Alkyl)		
Enzysept PLUS Wipes	Wool C&Tech	Quaternary Ammonium(N-Alkyl)		
Incidin OxyWipes	Ecolab	Hydrogen peroxide		
Metricide-OPA plus	Metrex	Ortho-phthalaldehyde		
Perasafe	LANXESS Deutschland GmbH	Citric acid		
Protex Ultra	Parker Laboratory Inc.	Quaternary Ammonium(N-Alkyl)		
SaFecide + (Sensitive wipes 200)	Liebecos	Quaternary Ammonium(N-Alkyl)		
Sani Cloth Active wipes	PDI	Quaternary Ammonium(N-Alkyl)		
Sani Cloth Universal wipes	PDI	Quaternary Ammonium(N-Alkyl)		
Sani-cloth bleach	PDI	Sodium hypochlorite	○	○
Sani-Cloth HB	PDI	Quaternary Ammonium(N-Alkyl)		
Sodium hypochlorite	Yuhan Chlorox	Hypochlorous acid		
Sono Ultrasound Wipes	Advanced Ultrasound Solutions Inc.	Quaternary Ammonium(N-Alkyl)	○	
Sporox II	DSHealthcare Inc.	Hydrogen peroxide		
Tristel duo	Tristel Solutions Limited	Chlorine dioxide		
Tristel Sporicidal Wipe	Tristel Solutions Limited	Chlorine dioxide		
T-spray II	Pharmaceutical Innovations, Inc.	Quaternary Ammonium(N-Alkyl)		
Virkon S	LANXESS	Pentapotassium, Succinic acid		
Wip Anios Excel	Laboratoires Anios	Quaternary Ammonium(N-Alkyl)		

※ Recommended by CDC (Centers for Disease Control and Prevention) for USA Market.  
CBC : Center for Biocide Chemistries.  
EPA : Environment Protection Agency


# Cautions & Warnings for disinfecting DR system

- Use only cleaning agents and disinfectants included in the list, especially for marked areas.
- Rub the product carefully with soft cloth when using disinfectants.



- Be careful not to let liquids or other substances get inside  
**: Do not spray directly** on the product, use a cloth or tissue.

- ※ Wear a face mask and gloves when cleaning or disinfecting the product.
- ※ Do not clean any electrical connectors, it can cause an electric shock or fire.



*Clean the surface in accordance with the disinfectant manufacturer's instructions (Temp, wiping/contact time etc.)*

# Disinfectant List for Detector All

· The following disinfectants can be used to disinfect detectors

: S-Detector™ (S4343-W, S4335-W, S3025-W), AccE Detector (S4343-AW, S4335-AW),  
AccE Standard Detector (S4343-AWM, S4335-AWM), AccE Glass-Free Detector (F4335-AW)

Name	Manufacture	Active Ingredient	CBC*	EPA*
Accel PREvention Wipes	Diversey Canada, Inc.	Hydrogen peroxide		
ALI-FLEX RTU	Lalema	Sodium hypochlorite		
ANIOSURF Premium NPC	ANIOS	Quaternary Ammonium(N-Alkyl)		
ANIOS社 surfa'safe	ANIOS	Quaternary Ammonium(N-Alkyl)		
Anioxyde 1000	Laboratoires Anios	Peracetic acid		
Azomax	Vernacare	5-chloro-2-methyl-2h-isothiazol-3-one		
Cleanisept Wipes	Dr.Schumacher GmbH	Quaternary Ammonium(N-Alkyl)		
Cleanisept wipe forte	Dr.Schumacher GmbH	Quaternary Ammonium(N-Alkyl)		
Cleanisept Wipes MAXI	Dr.Schumacher GmbH	Quaternary Ammonium(N-Alkyl)		
Clinell Sporidical wipes	GAMA Healthcare Ltd	Sodium Percarbonate		
Clinell Universal Wipes	GAMA Healthcare Ltd	Quaternary Ammonium(N-Alkyl)		
CLOROX Bleach Germicidal	The Chlorox Company	Sodium hypochlorite	○	○
CN'F	CARROLLASIA Co.,LTD	Alcohol		
Dettol	Reckitt Benckiser	Quaternary Ammonium(N-Alkyl)		
Distel Disinfectant Wipes	Tristel	Polymeric Biguanide Hydrochloride		
Ecolab Actichlor plus	Ecolab	Troclosene sodium		
ED Wipes	MH healthcare	Quaternary Ammonium(N-Alkyl)		
Enzysept PLUS Wipes	Wooil C&Tech	Quaternary Ammonium(N-Alkyl)		
Incidin OxyWipes	Ecolab	Hydrogen peroxide		
Lysoform	Lysoform	Formaldehyde		
Metricide-OPA plus	Metrex	Ortho-phthalaldehyde		
Mikrobac Tissues	BODE Chemie GmbH (HARTMANN)	Quaternary Ammonium(N-Alkyl)		
Perasafe	LANXESS Deutschland GmbH	Citric acid		
Protex Ultra	Parker Laboratory Inc.	Quaternary Ammonium(N-Alkyl)		
SaFecide + (Medical wipes 200)	Liebecos	Alcohol		
SaFecide + (Sensitive wipes 200)	Liebecos	Quaternary Ammonium(N-Alkyl)		
Sani-cloth bleach	PDI	Sodium hypochlorite	○	○
Sani-Cloth HB	PDI	Quaternary Ammonium(N-Alkyl)		
Sani-Cloth Super	PDI	Quaternary Ammonium(N-Alkyl)		
Sodium hypochlorite	Yuhan Chlorox	Hypochlorous acid		
Sono Ultrasound Wipes	Advanced Ultrasound Solutions Inc.	Quaternary Ammonium(N-Alkyl)	○	
Sporox II	DSHealthcare Inc.	Hydrogen peroxide		
Tristel duo	Tristel Solutions Limited	Chlorine dioxide		
Tristel Sporidical Wipe	Tristel Solutions Limited	Chlorine dioxide		
Trophon Companion wipe	Nanosonics Limited	Quaternary Ammonium(N-Alkyl)		
T-spray II	Pharmaceutical Innovations, Inc.	Quaternary Ammonium(N-Alkyl)		
Virkon S	LANXESS	Pentapotassium, Succinic acid		
Wip Anios Excel	Laboratoires Anios	Quaternary Ammonium(N-Alkyl)		

※ Recommended by CDC (Centers for Disease Control and Prevention) for USA Market.


CBC : Center for Biocide Chemistries.

EPA : Environment Protection Agency

# Cleaning Guide for S-Detector™ & AccE Detector


## Cleaning

- Always turn off the equipment completely before cleaning.
- Do not insert liquids into the system.
- When cleaning the exterior of the system, use a soft cleaning cloth dampened with tepid Water and soap.




## Liquid intrusion

- Do not spill liquid or chemicals onto the instrument or, in cases where the patient is injured, do not allow it to become wet with blood or other body fluids, as doing so may result in fire or electric shock.
- In such situation, protect the instrument with disposable covering as necessary.




S-Detector™ only




# Cleaning Guide for AccE Detector

- Recommended using a commercial plastic brush when cleaning non-slip pads.
- If need, non-slip pad can be exchanged or detached.



# Usage Guide for Detector Storage

- Detector has to be equipped properly.
- If you equip the detector in the storage rear for charging or moving with too much force or diagonally, pogo pin on the main cable connector may get damaged.
- There shouldn't be any moisture on the detector connector before equipping detector.
- If you equip the detector on the storage rear before completely removing moisture after cleansing, the moisture may get on the main cable connector, resulting in the damage of pogo pin.



# **Cleaning and Disinfecting the Ultrasound System System/Monitor/Probe**

**January, 2021  
Revision 10  
MI68-02691N\_H\_EN**

3366, Hanseo-ro, Nam-myeon, Hongcheon-gun, Gangwon-do, Republic of Korea  
[www.samsunghealthcare.com](http://www.samsunghealthcare.com)

Copyright © 2020 By SAMSUNG MEDISON CO., LTD. All rights reserved.



# Table of Contents

Disinfectants for System Surfaces (With Monitor) .....	4
Using Ultrasound Gel .....	8
Cleaning, Disinfecting, and Sterilizing the Probe.....	9
Validated High Level Disinfection Instruction for Samsung Transducers .....	12



## Disinfectants for System Surfaces (With Monitor)

The following disinfectant matrix indicates which disinfectants are compatible with specific system areas.

Please be cautious to avoid damage by only using approved disinfectants on areas listed.

No	Product Name	Active Ingredient	Manufacturer	Type	Body Surface	Control Panel	Monitor
1	Accel® INTERVention™ Wipes	Hydrogen peroxide	Diversey	Wipe	○	*	○
2	Accel® Prevention Wipes	Hydrogen peroxide	Diversey	Wipe	*	○	
3	ANIOSURF Premium NPC	Quaternary Ammonium (N-Alkyl)	Laboratoires Anios	Liquid	○	○	○
4	Avert Sporicidal disinfectant cleaner	Sodium hypochlorite	US Headquarters Diversey, Inc.	Liquid		○	
5	AZO detergent wipe	Isopropanol (70%)	Synergy Health Ltd	Wipe	*	○	
6	Azomax	2methyl-2h-isothiazol-3-one polyhexamethylene biguanide	Synergy Health Ltd	Wipe	*	○	
7	CLEANISEPT® WIPES FORTE	Quaternary Ammonium (N-Alkyl)	Dr.Schumacher GmbH	Wipe	○	*	○
8	Clinell Detergent Wipes (Yellow)	Lauryl Polyglucose, Butoxydiglycol	GAMA Healthcare Ltd	Wipe	○	*	
9	Clinell Sporicidal wipes	Citric acid	GAMA Healthcare Ltd	Wipe	○	○	
10	Clinell Universal Santitising wipes (Clinell Universal Wipes)	Quaternary Ammonium (N-Alkyl)	GAMA Healthcare Ltd	Wipe	*	○	○
11	Clorox Healthcare® Bleach Germicidal Wipes	Sodium Hypochlorite	The Clorox Company	Wipe	○	○	○
12	Clorox regular (Clorox regular1: water 10)	Sodium Hypochlorite	The Clorox Company	Liquid	○	○	○
13	CN'F	Alcohol	Neo Medical	Liquid	*	○	○

No	Product Name	Active Ingredient	Manufacturer	Type	Body Surface	Control Panel	Monitor
14	Dispatch@ Hospital Cleaner Disinfectant Towels with Bleach		The Clorox Company	Wipe	○	○	
15	Ecolab Actichlor Tablets	Adipic acid, Sodium Carbonate	Ecolab	Tablet	○	○	
16	ED Wipes	Quaternary Ammonium (N-Alkyl)	MH healthcare	Wipe	○	○	○
17	Enzysept PLUS Wipes	Quaternary Ammonium (N-Alkyl)	Wool C&Tech	Wipe	○	○	○
18	Ethanol 80%	Ethanol 80%	All manufacture	Liquid	*	○	
19	Incidin™ Oxy Wipe	Hydrogen peroxide	ECOLAB	Wipe	*	○	○
20	Meliseptol Wipes Sensitive	Quaternary Ammonium (N-Alkyl)	B. Braun Medical AG	Wipe	*	○	
21	Oxivir® Tb Wipes	Hydrogen peroxide	Diversey	Wipe	*	○	○
22	Perasafe	Citric acid	LANXESS Deutschland GmbH	Powder	○	○	○
23	PROTEX® ULTRA DISINFECTANT WIPES	Quaternary Ammonium (N-Alkyl)	Parker Laboratory Inc.	Wipe	○	○	○
24	saFecide Medical wipes	Alcohol	Liebecos	Wipe	*	○	○
25	saFecide Sensitive wipes	Quaternary Ammonium (N-Alkyl)	Liebecos	Wipe	*	○	○
26	Sani-Cloth® Active wipes		PDI	Wipe	○	*	
27	Sani-Cloth® bleach Germicidal Disposable Wipe	Sodium Hypochlorite	PDI	Wipe	○	○	○
28	Sani-Cloth® Detergent wipes		PDI	Wipe	*	○	
29	Sani-Cloth® HB Germicidal Disposable Wipe	Quaternary Ammonium (N-Alkyl)	PDI	Wipe	*	○	○

No	Product Name	Active Ingredient	Manufacturer	Type	Body Surface	Control Panel	Monitor
30	SciCan Optim 33TB Surface CleaneDisinfectant Wipes	Hydrogen peroxide	SciCan	Wipe	*	○	
31	Sono Ultrasound Wipes	Quaternary Ammonium (N-Alkyl)	Advanced Ultrasound Solutions Inc.	Wipe	*	○	○
32	Super Sani-Cloth® Germicidal Disposable Wipe	Quaternary Ammonium (N-Alkyl)	PDI	Wipe	*	○	○
33	Super Sani-Cloth® Plus Wipes		PDI	Wipe	*	○	
34	Surfa'safe	Quaternary Ammonium (N-Alkyl)	Laboratoires Anios	Liquid (Spray)	○	○	○
35	terralin® liquid	Alcohol	Schülke & Mayr GmbH	Liquid	*	○	
36	Tristel duo	Chlorine dioxide	Tristel Solutions Ltd.	Foam	*	○	○
37	Tristel Sporicidal Wipe	Chlorine dioxide	Tristel Solutions Ltd.	Wipe	*	○	○
38	Trophon® companion cleaning wipes	Quaternary Ammonium (N-Alkyl)	Nanosonics Limited	Wipe	*	○	
39	T-Spray II	Quaternary Ammonium (N-Alkyl)	Pharmaceutical Innovations, Inc.	Liquid (Spray)	○	○	○
40	Tuffie large detergent wipes	Ethylan CPG, Betain 810	Vernacare	Wipe	*	○	
41	Virkon S	Potassium peroxymonosulfate	Lanxess Corporation	Powder	○	○	
42	Wip'Anios Excel wipes	Quaternary Ammonium (N-Alkyl)	Laboratoires Anios	Wipe	○	○	○
43	Yuhanrox regular (Yuhanrox regular1: water10)	Sodium Hypochlorite	Yuhan clorox	Liquid	○	○	○

○	Compatible
*	If the disinfectant is used on non-approved area, it may result in surface cracking or discoloration.

**CAUTION:**

- ▶ Using an inappropriate disinfection procedure may damage the system. Make sure that you check the disinfectant's expiry date.
- ▶ Mix the disinfectant solution to the strength specified on the instruction of the disinfectant manufacturer.
- ▶ When using a disinfectant, wipe it carefully with a soft cloth.
- ▶ Air dry the device according to the disinfectant manufacturer's instruction, or remove water from all surface of the device using a sterile, lint-free cloth.
- ▶ Examine the device for damage such as cracks, splitting, sharp edges, or projections. If such damage is found, discontinue use of the device and contact your Samsung Medison representative.
- ▶ Repeated use of disinfectants may cause discoloration.

# Using Ultrasound Gel

For successful acoustic signal transmission, please only use the ultrasound gels approved by Samsung Medison.



## **WARNING:**

- ▶ The use of inappropriate ultrasound gels could result in damage to the probe. Using damaged probes may result in electric shocks and other hazards to the patients and/or users.
- ▶ Do not use ultrasound gels or coupling media that contain any of the following agents:
  - Oils such as mineral oil, cooking oil, gasoline, solvents, rust inhibitors, lanolin, paraffin-based grease, ester, and excessive silicon-based release agents;
  - Alcohols, such as acetone, methanol, and plasticizer (dioctylphtalate), or denatured alcohol;
  - Glacial acetic acid and iodine; or
  - All types of lotion or gel that contain aromatic substances



## **CAUTION:**

- ▶ When applying the ultrasonic gel to the probe, make sure that the tip of the gel tube does not touch the surface of the probe lens.
- ▶ The tip of the gel tube may cause damage to the probe lens.



# Cleaning, Disinfecting, and Sterilizing the Probe

All probes must be cleaned and disinfected after each use. Cleaning is an important procedure that must be carried out before disinfecting the probe. For information on cleaning and disinfecting the probe, please refer to 'Cleaning, Disinfecting, and Sterilizing the Probe' in the 'Probes' chapter of the user manual. Using an inappropriate disinfectant may damage the probe.



## WARNING:

- ▶ Always use protective equipment such as face mask, eyewear, and gloves when cleaning, disinfecting, and sterilizing probes.
- ▶ Inspect the housing, strain relief, lens and seal for damage, and check for any functional degradation after cleaning and disinfecting the probe.
- ▶ Using an inappropriate cleaning or disinfecting agent may damage the probe.

## Information on Detergent, Disinfectant, and Ultrasound Gel

### Reprocessing Method by Probe Type

To maintain the performance of ultrasound probes, proper maintenance is required.

As ultrasound probes are classified into critical, semi-critical or non-critical devices based on the standards of FDA guidance\* and the Hygiene Requirements for the Reprocessing of Medical Devices from Germany guideline of Robert Koch Institute (RKI), proper cleaning and disinfection or sterilization methods for that classification should be used. Therefore, you should use the cleaning, disinfection, and sterilization methods appropriate for each classification.

#### Choosing the Correct Probe-Care Method in below Table

Classification Criteria	Contact Area	Application Probe	Level Selection
Non-critical medical devices	Intact skin	Curved, Linear, and Phased array probes	Low level disinfection
Semi-critical medical devices	Mucous membrane, damaged skin	Endocavity, MPTEE	High level disinfection or sterilization
Critical medical devices	Blood, sterile tissue, etc.	Intraoperative	Sterilization

\* Guidance for Industry and FDA Staff – Marketing Clearance of Diagnostic Ultrasound Systems and Transducers - Appendix E

\* The FDA reprocessing guidance 'Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling, Guidance for Industry and Food and Drug Administration Staff' March 17, 2015, (<https://www.fda.gov/media/80265/download>)

The care method for your probes determines the appropriate disinfectant for your probe. An appropriate detergent, disinfectant or ultrasound gel should be used for all probes. For details about compatible detergent, disinfectants, and ultrasound gel, please see 'Disinfectants Matrix' on Samsung Medison website and the User Guide.

- ▶ User Guide: This is provided as a booklet upon purchase of the product.







# Validated High Level Disinfection Instruction for Samsung Transducers

## ✚ Tristel Duo ULT: Wiping Method for Endocavity Probes

### 1. Cleaning the Probe




**NOTE:** Always use protective equipment such as face mask, eyewear, and gloves when cleaning and disinfecting probes.

- 1) If applicable, remove the biopsy adapter, needle guide or protective sheath from the probe.
- 2) Before conducting the HLD process, clean the probes thoroughly.



**NOTE:** Do not use paper or abrasive products when cleaning the probe. They can damage the lens of the probe.

- 3) Hold the probe and put one pump of Tristel Duo ULT in the crevices of the probe. (See Figure 1 ❶)
- 4) Use a swab to spread the Tristel Duo ULT Foam in the crevices. (See Figure 1 ❷)
- 5) Put two pumps of Tristel Duo ULT on a Tristel Dry Wipe and clean the transducer from housing to the lens of the probe. (See Figure 1 ❸)
- 6) Use another dry swab to dry and clean the crevices. (See Figure 1 ❹)
- 7) After cleaning, thoroughly dry the probe.




[Figure 1]

### 2. Using Tristel Duo ULT


Follow Tristel Duo ULT manufacturer's instructions for preparation, temperature, solution strength, and duration of contact. Ensure that the solution strength and duration of contact are appropriate for the intended clinical use of the device. Be sure to observe the solution expiration date.

### 3. Disinfecting the probe

- 1) Change your gloves.
- 2) Hold the probe and put one pump of Tristel Duo ULT in the crevices of the probe.  
(See Figure 2 ❶)
- 3) Use a swab to spread the Tristel Duo ULT Foam in the crevices. Leave the disinfectant on the probe for at least 30 seconds. (See Figure 2 ❷)
- 4) Put two pumps of Tristel Duo ULT on a Tristel Dry Wipe and disinfect the transducer from housing to the lens of the probe. Leave the disinfectant on the probe for at least 30 seconds. (See Figure 2 ❸)



[Figure 2]



[Figure 3 – Crevices]

### 4. Dry

- 1) Leave the surface to dry to ensure at least 30 second contact time.
- 2) Use the probe immediately or store the probe properly (e.g. in a sterile bag) to prevent any further contamination.

### 5. Management

- 1) Register the disinfection data for traceability according to local regulations.
- 2) Examine the probe and cable for damage such as cracks, splitting, sharp edges, or projections. If such damage is found, discontinue use of the device and contact your Samsung Medison representative.

※ Using disinfection wipes for devices requiring high-level disinfection may be prohibited in some countries such as USA.

## ✚ Metricide OPA Plus: Immersion Method for Endocavity Probes

### 1. Cleaning the Probe



**NOTE:** Always use protective equipment such as face mask, eyewear, and gloves when cleaning and disinfecting probes.

- 1) If applicable, remove the biopsy adapter, needle guide or protective sheath from the probe.
- 2) Before conducting the HLD process, clean the probes thoroughly with wiping or immersion method as per the instructions of the chosen detergent manufacturer.



**NOTE:** Do not use paper or abrasive products when cleaning the probe. They can damage the lens of the probe.

- 3) If necessary, follow the rinse (and neutralization) instructions from the cleaner or disinfectant manufacturer for rinsing the probe.
- 4) After cleaning, thoroughly dry the probe.

### 2. Using Metricide OPA Plus

Follow Metricide OPA Plus manufacturer's instructions for preparation, temperature, solution strength, and duration of contact. Ensure that the solution strength and duration of contact are appropriate for the intended clinical use of the device. Be sure to observe the solution expiration date.


### 3. Disinfecting the Probe

- 1) Change your gloves.
- 2) Immerse the probe into the Metricide OPA Plus solution. (See Figure 4)



**NOTE:** Ensure that all air bubbles are removed from the surface of the probe by agitating.

- 3) Allow the probe to soak in the Metricide OPA solution for 12 minutes.



[Figure 4]

#### **4. Rinse**

- 1) Thoroughly rinse the probe by fully immersing it in pure water to the point in the figure above. (See Figure 4)
- 2) Remove all air bubbles from the surface of the probe and allow the probe to set for a minimum of 1 minute.
- 3) Repeat the Step 1) and 2) of Rinse two more times for a total of 3 rinses using a fresh batch of pure water each time.

#### **5. Dry**

- 1) Remove water from all surfaces of probe with a sterile lint-free cloth.
- 2) Use the probe immediately or store the probe properly (e.g. in a sterile bag) to prevent any further contamination.

#### **6. Management**

- 1) Register the disinfection data for traceability according to local regulations.
- 2) Examine the probe and cable for damage such as cracks, splitting, sharp edges, or projections. If such damage is found, discontinue use of the device and contact your Samsung Medison representative.

## ❏ Trophon EPR, Trophon 2: Disinfection System Method for Endocavity Probes

### 1. Cleaning the Probe



**NOTE:** Always use protective equipment such as face mask, eyewear, and gloves when cleaning and disinfecting probes.

- 1) If applicable, remove the biopsy adapter, needle guide or protective sheath from the probe.
- 2) Before conducting the HLD process, clean the probes thoroughly with wiping or immersion method as per the instructions of the chosen detergent manufacturer.



**NOTE:** Do not use paper or abrasive products when cleaning the probe. They can damage the lens of the probe.

- 3) If necessary, follow the rinse (and neutralization) instructions from the cleaner or disinfectant manufacturer for rinsing the probe.
- 4) After cleaning, thoroughly dry the probe.

### 2. Preparing Trophon EPR, Trophon 2




**NOTE:** Refer to the user instructions of the Trophon EPR, Trophon 2 for details of proper use, storage and disposal of waste.


- 1) Place the probe in Trophon EPR, Trophon 2 following the user instructions of Trophon EPR, Trophon 2.



**NOTE:** Improper positioning of the probe can damage the probe due to the increased temperature of the wall and high-level disinfection may not be achieved. (See Figures 5, 6)



[Figure 5 – Trophon EPR]



[Figure 6 – Trophon 2]

### 3. Disinfecting the Probe

- 1) Change your gloves.
- 2) Start the disinfection process and wait until the end of the cycle. (The standard cycle takes 7 minutes)
- 3) Wear a new set of clean gloves and always remove the probe immediately at the end of the cycle from the Trophon EPR, Trophon 2.



**NOTE:** After completion of the disinfection process, the surface temperature of the probe and chamber may have increased. Always use protective equipment such as face mask, eyewear, and gloves.

### 4. Rinse & Dry

- 1) Remove hydrogen peroxide from all surfaces of the probe with a sterile lint-free cloth.



**NOTE:**

- ▶ Aggressive wiping or scrubbing can damage the probe. Use a gentle wiping motion for the probe lens, the strain relief, and the areas surrounding the strain relief.
- ▶ When wiping the probe, hold the housing. Do not suspend the probe by the cable. This can damage probes.

- 2) Use the probe immediately or store the probe properly (e.g. in a sterile bag) to prevent any further contamination.

### 5. Management

- 1) Register the disinfection data for traceability according to local regulations.
- 2) Examine the probe and cable for damage such as cracks, splitting, sharp edges, or projections. If such damage is found, discontinue use of the device and contact your Samsung Medison representative.

## Cidex OPA: Immersion Method for MPTEE Probe

### 1. Cleaning the Probe



**NOTE:** Always use protective equipment such as face mask, eyewear, and gloves when cleaning and disinfecting probes.

- 1) If applicable, remove the biopsy adapter, needle guide or protective sheath from the probe.
- 2) Before conducting the HLD process, clean the probes thoroughly with wiping or immersion method as per the instructions of the chosen detergent manufacturer.



**NOTE:**

- ▶ Do not use paper or abrasive products when cleaning the probe. They can damage the lens of the probe.
- ▶ Immerse only the portion of the probe between the tip at the end of the probe and the area marked as 100cm. (See Figure 7)

- 3) If necessary, follow the rinse (and neutralization) instructions from the cleaner or disinfectant manufacturer for rinsing the probe.
- 4) After cleaning, thoroughly dry the probe.

### 2. Using Cidex OPA

Follow Cidex OPA manufacturer's instructions for preparation, temperature, solution strength, and duration of contact. Ensure that the solution strength and duration of contact are appropriate for the intended clinical use of the device. Be sure to observe the solution expiration date.


### 3. Disinfecting the Probe

- 1) Change your gloves.
- 2) Immerse the probe into the Cidex OPA solution. (See Figure 7)



**NOTE:** Immerse only the portion of the probe between the tip at the end of the probe and the area marked as 100cm.





[Figure 7]



**NOTE:**

- ▶ Ensure that all air bubbles are removed from the surface of the probe by agitating.
- ▶ Do not allow sharp objects to touch the lens. Do not touch the lens unnecessarily. This can damage the lens.

- 3) Allow the probe to soak in the Cidex OPA solution for 12 minutes.

#### 4. Rinse

- 1) Thoroughly rinse the probe by fully immersing it in sterile pure water to the point in the figure above. Immerse only the portion of the probe between the tip at the end of the probe and the area marked as 100cm. (See Figure 7)
- 2) Remove all air bubbles from the surface of the probe and allow the probe to set for a minimum of 1 minute.
- 3) Repeat the Step 1) and 2) of **Rinse** two more times for a total of 3 rinses using a fresh batch of sterile pure water each time.

#### 5. Dry

- 1) Remove water from all surfaces of probe with a sterile lint-free cloth.
- 2) Use the probe immediately or store the probe properly (e.g. in a sterile bag) to prevent any further contamination.

#### 6. Management

- 1) Register the disinfection data for traceability according to local regulations.
- 2) Examine the probe and cable for damage such as cracks, splitting, sharp edges, or projections. If such damage is found, discontinue use of the device and contact your Samsung Medison representative.

## ✚ Gigasept® PAA concentrate: Immersion Method for MPTEE Probe

### 1. Cleaning the Probe



**NOTE:** Always use protective equipment such as face mask, eyewear, and gloves when cleaning and disinfecting probes.

- 1) If applicable, remove the biopsy adapter, needle guide or protective sheath from the probe.
- 2) Before conducting the HLD process, clean the probes thoroughly.



**NOTE:**

- ▶ Do not use paper or abrasive products when cleaning the probe. They can damage the lens of the probe.
- ▶ Immerse only the portion of the probe between the tip at the end of the probe and the area marked as 100cm. (See Figure 8)

- 3) Wipe the MPTEE probe with 5 Tristel Dry Wipes soaked with 8 ml 1.6% CIDEZYME and wipe for 3 minutes.
- 4) Thoroughly clean the probe with a detergent especially taking care of the junctions of the flexible part and the crevices on the transducer. (See Figure 8)
- 5) Rinse the MPTEE probe for 1 minute with water.



[Figure 8]

### 2. Using Gigasept® PAA concentrate

Prepare a Gigasept® PAA concentrate with a concentration of at least 2% using the manufacturer's instructions.

### 3. Disinfecting the Probe

- 1) Change your gloves.
- 2) Place the MPTEE probe in the 2% Gigasept® PAA concentrate. (See Figure 9)




**NOTE:** Immerse only the portion of the probe between the tip at the end of the probe and the area marked as 100cm. (See Figure 9)

- 3) Ensure that all air is removed from the surface of the MPTEE probe.
- 4) Leave the MPTEE probe for 15 minutes in the solution.



**NOTE:** Gigasept® PAA concentrate shall not be reused beyond 12 hours.



[Figure 9]

#### 4. Rinse

- 1) Remove the MPTEE probe from the solution and thoroughly rinse all surfaces up to the immersion point. (See Figure 9)
- 2) Repeat this process two more times with sterile water, each time for a minimum of 1 minute.

#### 5. Dry

- 1) Thoroughly dry all surfaces of the probe using a sterile, lint-free wipe or cloth, changing wipes/cloths when necessary to ensure the probe is completely dry. Visually inspect the probe to ensure all surfaces are clean and dry. Repeat the drying steps if any moisture is still visible.
- 2) Use the probe immediately or store the probe properly (e.g. in a sterile bag) to prevent any further contamination.

#### 6. Management

- 1) Register the disinfection data for traceability according to local regulations.
- 2) Examine the probe and cable for damage such as cracks, splitting, sharp edges, or projections. If such damage is found, discontinue use of the device and contact your Samsung Medison representative.